




ON SITE DRUG SCREENING DEVICES

KahnTact USA, Inc. On-Site Drug Screens utilize the same technology that laboratories use. Similar to the laboratories initial screening process the tests utilize enzyme immune assays to detect the use of illicit drugs in a urine specimen. These On-Site Drug Screen Devices provide results in a few minutes without the need for expensive machinery or technical training.

All of our tests are easy to use and read. Just add the test to urine specimen and read the results. The appearance of one test result line indicates a positive result. Two lines

indicate a negative result. The tests are FDA cleared and accurately set to NIDA SAMHSA cutoff levels.



Tests can be configured in a variety of combinations. Individual drug screens are available in dip/read or drop/read formats. There are 10 different drugs to choose from to configure a multi-panel test to suit your company's needs:

AVAILABLE DRUG METABOLITES	DETECTION LEVEL
(AMP) Amphetamine.....	1000 ng/mL
(MET) Methamphetamine.....	1000 ng/mL
(BZO) Benzodiazapines.....	300 ng/mL
(MTD) Methadone.....	300 ng/mL
(BAR) Barbiturates.....	300 ng/mL
(THC) Marijuana.....	50 ng/mL
(PCP) Phencyclidine.....	25 ng/mL
(OPI) Opiate.....	2000 ng/mL
(MOP) Opiate.....	300 ng/mL
(COC) Cocaine.....	300 ng/mL
(TCA) Tricyclic.....	1000 ng/mL

INTEGRATED CUP



The integrated cups work in the same way as other on-site drug screens. The cups, however, are fully integrated. This alleviates the possibility of the test administrator coming into contact with the urine specimen. Furthermore, the sample can be separated to send off to a lab for confirmation of positive results. The lock mechanism eliminates any contamination of the sample.

The small size reduces shipping costs and decreases storage requirements. The cup has an easy Xerox feature where the cup can be laid down on its side and copied for retention of results.

ALCOHOL SALIVA TEST

The alcohol screen is a one-step rapid test that utilizes a highly sensitive method to detect alcohol presence. The tests are available in a semi-quantitative or qualitative approximation of blood alcohol levels.

- * Rapid 2 minute test
- * Color gradient indicates .02-.3 alcohol levels
- * Available in boxes of 24
- * D.O.T. Compliant



LAB BASED SUBSTANCE ABUSE TESTING

KahnTact USA, Inc. provides third party administration services for several of the largest clinical laboratory networks in the United States including LabCorp and Quest Diagnostics. Utilizing several labs has allowed us to construct one of the most comprehensive clinical test menus in the occupational and clinical testing industry. The combination of several labs gives our clients effortless access to thousands of Patient Service Centers nationally.

KahnTact USA, Inc. provides routine occupational and clinical specimen collections for testing and screening. In addition to routine urine drug screens, we offer procedures such as blood specimen collections, paternity collections and breath alcohol testing. Equally important types of drug testing can be done including, random testing, voluntary testing, reasonable suspicion testing, specific condition testing and follow-up testing.

DRUG TEST PANELS

In addition to the standard NIDA (National Institute Drug Administration) panel, KahnTact USA, Inc. offers a variety of test panels and cut off levels, which can be combined with adulterant testing and MRO (Medical Review). Approximately 95 percent of negative results are reported within 24 hours of receipt into our laboratories. The remaining 5 percent of negatives and positives are reported within 72 hours of receipt by our laboratories. A list of the most commonly drugs tested:

Amphetamines	Cannabinoids
Cocaine	Opiates metabolites
Phencyclidine	Benzodiazapines
Barbiturates	Propoxyphene
Methadone	Methaqualone
Alcohol	Ecstasy

ADULTERANT TESTING

Adulterant testing detects substances (either added to the urine specimen or ingested) that are promoted as “cleansing agents” to prevent the detection of drug use. Adulterant tests are designed to detect specimen tampering or substitution and provide the best information from which to make informed personnel decisions. The adulterants tested for include these compounds:

Nitrite > 500 ug/mL	Pyridinchlorochromate
Creatinine < 20 mg/dL	Acids, Bases
PH(<) or (>11)	Bleach
Specific Gravity indicated	Glutaraldehyde
Alternate Screening indicated	Chromates



In addition to pre-employment and random drug screens, KahnTact USA, Inc. is prepared to be your one-source provider for a variety of ancillary services, including physical examination management, workplace exposure monitoring, other specimen collection, and scheduling services.

NATIONAL CALL CENTER

KahnTact USA, Inc. also provides a skilled staff who can respond immediately to your inquiries. KahnTact USA, Inc. is dedicated to resolving inquiries precisely and completely.

COLLECTIONS MANAGEMENT

Through its partner networks KahnTact USA, Inc. offers drug screen collections at more than 2000 patient service centers located in the United States. The network of approved sites provides services including, urine drug screen collections, breath alcohol testing, on-site testing products collections, physical examination and on-location collections. Information on site locations is accessible through our toll free number or online.

**POST ACCIDENT COLLECTION
COORDINATION PROGRAM**

When an accident happens after hours, a laboratory representative is only a phone call away to assist in locating a walk-in site or mobile collector, 24 hours a day, 7 days a week, 365 days a year.

REPORTING

KahnTact USA, Inc. offers multiple results-delivery systems to meet specific customer needs. Available reporting options for lab results include:

- Automatic faxing
- Voice reporting (Spanish available)
- On-line specimen status inquiry